

Module 3: Immunostaining

BMES Cell Team

Fall 2020



Outline

- Review of Modules 1 and 2
- Introduction to Immunostaining
- Antibody-Antigen Reaction
- Overview of the Protocol

Serial Dilutions for Western Blotting

- When you're dealing with changes in **concentration**, the two most important equations are:

1. Beer's Law

$$A = \epsilon \ell C \quad (1)$$


2. The Dilution Equation

$$C_1 V_1 = C_2 V_2 \quad (2)$$

Beer's Law


$$A = \epsilon \ell C$$

- A is the Absorbance
- ϵ is a constant that depends on the substance we are dealing with
- ℓ is the length of the container or path
- C is the concentration of the substance we are dealing with




Beer's Law

- If the solution is concentrated, more light would be absorbed.




- If the path length ℓ is long, more light would also be absorbed.



Serial Dilutions for Western Blotting

Demonstration



Lab Protocol

Part 1: Serial Dilutions

1. Pipet 20 μL of H_2O from B1 down to G1
2. Pipet 20 μL protein standard into B1, mix well, then bring 20 μL down into the next and mix well
3. Continue performing serial dilutions until the last well

Part 2: Creating Triplicates

1. Starting from the least concentrated standard in G1, transfer 5 μL to each of the three wells on the right
2. Continue upwards until all replicates have been made


Interpreting Western Blot Data

- In Module 2, we derived one of the kinematic equations to show the relationship between mass and displacement:

$$x = \frac{F}{2m} t^2 + v_0 t + x_0 \quad (3)$$


- From this equation, we see that the mass of a protein (m) varies inversely with the displacement
- Thus, smaller proteins will travel further in Gel Electrophoresis

Interpreting Western Blot Data



Modelling Cell Growth with Differential Equations

- Don't worry about knowing how to solve the differential equation since I just wanted you to know where it is mathematically derived from
- However, the point I want you to take away from this is that most microorganisms (i.e. cells, bacteria, viruses) exhibit a **logistic growth pattern**, which can be modelled thru differential equations



Introduction to Immunostaining


- **Definition:** **Immunostaining** is an antibody-based method used to detect the presence of a specific protein in cells and tissues.
- Note that this is a **qualitative** procedure, since we are only trying to detect the presence of a specific protein
 - This is different from Western Blotting because we don't care about the exact concentration of a protein

Applications of Immunostaining

- **Cancer Detection** We can use immunostaining to detect prognostic markers for cancer.
- **Viral Infections** We can use antibodies specific to a virus to test whether a patient is indeed infected.
- **Genetic Testing** If a genetic disease is widespread in your family, you can use immunostaining to determine whether you will develop that disease.

Antibody-Antigen Reaction

Parts of the Antibody



Overview of the Protocol

1. Fix cells
2. Permeabilize the cell membrane
3. Incubate with primary antibodies
4. Incubate with secondary antibodies
5. Image cells with microscope

The Experiment

- If campus reopens this year, we will perform this experiment using the following cells:
 - **3T3 Mouse Fibroblasts**
 - **Neural Stem Cells**
- Also, we will use immunostaining to detect the following proteins:
 - **Actin** — A protein found in almost all cells and makes up the cytoskeleton. It is used to drive movement and cellular transport processes.
 - **SOX2** — A transcription factor that helps maintain the *pluripotency* of stem cells.

1. Fixing Cells

What is the purpose of fixing cells?


- It halts all cellular processes and immobilizes proteins
- This makes it easier for the antibodies to access and bind to proteins

How do we fix cells?


- We use a 4% solution of paraformaldehyde (PFA)
 - This reagent allows proteins to form covalent bonds with each other and their surroundings

1. Fixing Cells


Unfixed Cell



Fixed Cell



1. Fixing Cells





Think of it as popping a balloon filled with **water** vs. a balloon filled with **ice**.

- To make it easier to detect antibodies, we want the cell to stay intact after we lyse the membrane
 - Thus, we prefer the ice-filled balloon

1. Fixing Cells

Here is a layout of the 24-well plate that we will be using:




2. Permeabilization

- Intercellular vs. Membrane-Bound Proteins:
 - Based on their functions described earlier, which category do you think Actin and SOX2 fit under?
- For *intercellular* proteins, we need to lyse the cell membrane so the antibodies can access them
- In order to accomplish this, we use a reagent called **PBS-Triton**
 - PBS-Triton removes the membrane-bound proteins, which are responsible for keeping the membrane intact
- Also, we need to wash with **PBS-Tween 20** before and after lysis
 - PBS-Tween 20 removes background noise caused by non-specific binding, thus helping you obtain more accurate results


2. Permeabilization

- Let's take a closer look at PBS-Triton:



2. Permeabilization

Protein-Surfactant Complex




**Surfactant
(PBS-Triton)**



3. Primary Antibody Incubation

- The primary antibody we will use is **Rabbit anti-SOX2**
 - This should only bind to NSCs and not 3T3s
- Incubate overnight



3. Primary Antibody Incubation

- We need to add 500 μL of antibody solution to each well, and we have a total of ten wells. What is the total volume of solution we need?

$$(500 \mu\text{L}) \times (10 \text{ wells}) = 5 \text{ mL}$$

- The SOX2 solution should be diluted at a factor of 1:200. To prepare this solution with 5 mL of PBS-Tween 20, how much SOX2 should we add?

$$25 \mu\text{L SOX2}$$

4. Secondary Antibody Incubation

- Again, we want to block with PBS-Tween 20 to prevent nonspecific binding
- All secondary reagents will be wrapped in foil to avoid photobleaching
 - If reagents are exposed to light, they won't be able to fluoresce when we image it later
- The secondary antibodies used will be:
 - **Hoechst** — Used as a control to stain all nuclei to fluoresce blue
 - **Goat anti-Rabbit** — Binds to Rabbit anti-SOX2 to fluoresce green
 - **Phalloidin** — Binds to Actin in the cytoskeleton to fluoresce red
- We want to incubate for one hour

4. Secondary Antibody Incubation

- Again, we need to add 500 μL of antibody solution to each well, and we have a total of ten wells.
 - Thus, we need 5 mL of total solution
- The Hoechst and Goat anti-Rabbit solution should both be diluted at 1:1000. Phalloidin should be diluted at 1:500. How much antibody is required to prepare each solution if you start out with 5 mL of PBS-Tween 20?

5 μL Hoechst
5 μL Goat anti-Rabbit
10 μL Phalloidin

5. Imaging

- Again, turn off the lights to avoid photobleaching
- Use a transmitted channel at 4x and 10x objective to view the cells
- Use a fluorescent channel to view the secondary antibodies:
 - For **Hoechst**, use DAPI at 10x
 - For **SOX2**, use GFP at 10x
 - For **Phalloidin**, use RFP at 40x

