


### Protocol Discussion 3: Cell Viability Analysis


In Tuesday's protocol walkthrough, we discussed a cell viability experiment where 3T3 cells were incubated with different concentrations of DMSO and stained with calcein and ethidium homodimer (Live-Dead Stain). In this protocol discussion, we will be utilizing ImageJ to quantify the number of live and dead cells and differences in their morphology.

The figure below depicts 3T3 cells treated with 2%, 5%, 10%, and 25% DMSO under the GFP and RFP fluorescent microscope filters.


2% DMSO




5% DMSO



10% DMSO



25% DMSO



## Overview of Manual Cell Counting

1. Download the .jpg files from the Cell Team site
2. Open the “2% DMSO Treatment” .jpg file in ImageJ
3. Resize the image to the area that you want to work with
4. Calibrate the pixel size to the 200 um bar in the image using the “Set Scale” feature
5. Use the multipoint tool to mark live cells in yellow and record the live cell number in the table
6. To count the dead cells, you can change the multipoint counter color, re-open the image and use the multipoint counter, or just count them without the multipoint tool (recommended method for this set of data, as there are so few dead cells)
7. Repeat for the other three images (you can divide up this work in breakout rooms)

	Number of Live Cells	Number of Dead Cells	% Live Cells
<b>2% DMSO Treatment</b>			
<b>5% DMSO Treatment</b>			
<b>10% DMSO Treatment</b>			
<b>25% DMSO Treatment</b>			

**How did increasing concentration of DMSO impact cell viability?**

---

---

---

**Researchers have studied cancer cells in low concentrations (<1%) of DMSO to study its effect as a chemotherapeutic. Why do you think researchers have been exploring the use of DMSO? Based on your finding above, what are some of your concerns with using DMSO as a chemotherapeutic?**

---

---

---

---

---

---

## Quantifying Cell Morphology

1. Open the “2% DMSO Treatment” .jpg file in ImageJ
2. Select the area and perimeter measurements
3. Zoom into the image using the magnifying glass tool
4. Select the freehand tool and trace one of the live cells
5. Click Analyze -> Measure to record the area and perimeter of that cell
6. Click Process -> Noise -> Add Noise. This will help you keep track of the cells you have already counted
7. Repeat Steps 5-6 for all the live cells in the image
8. Copy-paste the live cell data into an Excel Sheet
9. Compute the circularity of all the cells, average circularity, average area, and average cell perimeter, and record your values below.
  - Note:  $\text{Circularity} = 4 \cdot \pi \cdot (\text{area} / \text{perimeter}^2)$
  - A circularity of one means that the cell is a perfect circle
10. Repeat Steps 5-9 for all the dead cells in the image
  - Note: this process is very tedious, so divide up the work in your breakout rooms!

LIVE CELLS	Average Area	Average Perimeter	Average Circularity
2% DMSO Treatment			
5% DMSO Treatment			
10% DMSO Treatment			
25% DMSO Treatment			

DEAD CELLS	Average Area	Average Perimeter	Average Circularity
2% DMSO Treatment			
5% DMSO Treatment			
10% DMSO Treatment			
25% DMSO Treatment			

**Compare your live cell and dead cell average area, perimeter, and circularity measurements. Why do you think the circularities are different? (Note: more circular cells have a circularity measurement closer to 1, and less circular cells have a circularity measurement closer to zero).**

---

---

---

---

---

---

---

---

---

---

**What are the limitations of cell counting for this particular sample using this method? What modifications can be made to the procedure to enhance the data acquisition process and quality of the data?**

---

---

---

---

---

---